

DIPYRIDYL Paper

for the rapid determination of ferrous ions

Colour reaction:

In the presence of Fe(II)-ions a red spot appears on white background.

Method of application:

Apply a drop of the weak mineral acid test solution (> pH 2) to the test paper. In the presence of Fe(II)-ions a red spot appears, or in the presence of small quantities of Fe(II)-ions a red ring.

Limit of sensitivity: 2 mg/L Fe(II)

Interferences:

The paper is specific for the determination of Fe(II)-ions. A few other metals result in pale discolourations of the test paper, but these are covered by the colour reaction of the ferrous ions. Fe(III)-ions result in a pale yellow reaction. In the presence of substantial quantities of Fe(III)-ions and small quantities of Fe(II)-ions, the Fe(III)-ion can be masked by the addition of small quantities of NaF.

CTL SCIENTIFIC SUPPLY CORP.

1016-3 Grand Boulevard, Deer Park, NY 11729

Tel: **631-242-4249**

Web: **www.ctlscientific.com**

Manufacturer: **Macherey-Nagel GmbH & CO. KG**

Rev: 2025-03