

COBALT Test Paper

for the rapid determination of cobalt

Colour reaction:

The white test paper turns blue upon contact with cobalt(II)ions.

Method of application:

Apply a drop of the mineral acid solution ($\text{pH}<1$) to be tested to the test paper. The presence of substantial quantities of cobalt(II)ions is indicated by the appearance of a blue spot; smaller amounts result in a blue ring.

Limit of sensitivity: 25 mg/l Co^{2+}

Interferences:

Iron and copper interfere.

Ferric-ions in substantial quantities result in the appearance of a reddish-brown spot, which, however, disappears again after a few seconds.

Cupric-ions show up as a brown ring. After 3 to 5 minutes, small amounts of cobalt combined with substantial quantities of copper are detected through the appearance of the blue cobalt ring within the brown copper ring, partly overlapping the latter. In doubtful cases, a drop of pure cupric-solution can be applied next to the test solution.

Product data and ordering information

Type	Qualitative test paper
Presentation	Pack of 100 strips 20x70 mm
Color reaction	white to blue
Gradation	25 mg/l (ppm) Co^{2+}
REF	90728

CTL SCIENTIFIC SUPPLY CORP. 1016-3 Grand Boulevard Deer Park, NY 11729
Tel: 631-242-4249

web: www.ctlscientific.com

Manufacturer: Macherey-Nagel GmbH & CO. KG