

BISMUTH Test Paper

for the rapid determination of bismuth

Colour reaction:

In the presence of Bi^{3+} , the test paper shows orange-red spots against a pale-yellow background.

Presentation:

Plastic boxes of 200 strips measuring 20 x 70 mm each.

Method of application:

Apply a drop of the hydrochloric acid solution to the test paper. In the presence of Bi^{3+} an orange-red spot appears on the paper, in the case of small quantities of Bi^{3+} an orange-red ring. Since, in the presence of minute quantities of Bi^{3+} , the color reaction is not permanent, the determination is to be made immediately upon application of the test solution to the paper.

Yellow spots do not indicate a positive reaction.

Limit of sensitivity: 60 mg/l Bi^{3+}

Interferences:

A number of metal ions produce a yellow color reaction, which, however, does not interfere with the determination of bismuth.

Ferric ions in substantial amounts produce a clear, black ring next to the reaction spot for Bi^{3+} .

Nitric acid limits the sensitivity of the test paper. The test solution should not contain more than 2-3% HNO_3 .

Product data and ordering information

Type	Qualitative test paper
Presentation	Pack of 200 strips 20x70 mm
Color reaction	pale-yellow to orange-red
Gradation	60 mg/l Bi^{3+} (Nitric acid >3% leads to reduced sensitivity)
REF	90733

CTL SCIENTIFIC SUPPLY CORP. 1016-3 Grand Boulevard Deer Park, NY 11729

Tel: 631-242-4249

web: www.ctlscientific.com

Manufacturer: Macherey-Nagel GmbH & CO. KG