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Instructions for use Indicator paper for udder control

This indicator paper is used for the rapid detection of a disturbed secretion, which can often be interpreted as a first sign of a beginning udder disease.

The paper contains four spots – corresponding to the four quarters of the udder – which are impregnated with a special indicator. From each teat a thin stream of the initial milking is applied to the corresponding indicator spot.

For healthy udders, all spots assume the same yellowish-green colour, while a disease is indicated by deviating colours of the individual quarters. Severely sick quarters are indicated by a bluish-green to blue colour, while slightly infected quarters give green spots. Milk from fresh-milking cows can sometimes result in a quince-yellow colour, even for healthy cows. The main focus of the evaluation should be directed on a comparison of the colours of the four spots i.e. whether they are equal, or whether an individual spot deviates from the others.


The air in a milking parlour contains ammonia and may change the colouring of the test paper to light green – therefore keep the paper in a closed airtight bag outside the milking parlour.

CTL SCIENTIFIC SUPPLY CORP. 1016-3 Grand Blvd. Deer Park, NY 11729

Tel: 631-242-4249

web: www.ctlscientific.com

Manufacturer: Macherer-Nagel GmbH & CO. KG

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