

QUANTOFIX® Glucose

Pack content:

1 aluminum container with 100 test strips

Measuring range:

Visually	Reflectometrically
50–2000 mg/L or 0.3–11.1 mmol/L	50–2000 mg/L or 0.3–11.1 mmol/L

Color gradation:

0 · 50 · 100 · 250 · 500 · 1000 · 2000 mg/L or

0 · 0.3 · 0.55 · 1.4 · 2.8 · 5.55 · 11,1 mmol/L

Hazard warnings:

This test does not contain hazardous substances that must be labelled.

General indications:

Remove only as many test strips as are required. Close the container immediately after removing a strip.

Do not touch the test field.

Instructions for use:

1. Dip test strip into test solution for 2 s.
 2. Shake off excess liquid.
 3. Wait 30 s.
 4. Compare with the color scale. If Glucose is present, the test field turns green-blue.
- The reaction color of the test field may change after the value has been taken. It is therefore crucial to evaluate the coloration within the prescribed time scale in order to achieve a correct result. Ignore color changes that occur after the reaction time (30 seconds).

Interferences:

If the sample solution has a pH value of 2–9, the reaction will take place without interferences. Strongly acidic solutions can be adjusted with sodium hydroxide solution, alkaline solutions with hydrogen chloride solution to a pH of 5–7. The following ions interfere with the determination when the concentrations below are exceeded (tested with 0 and 250 mg/L respectively of glucose solutions):

1000 mg/L: Ca^{2+} , Cl^- , HCO_3^- , NO_3^- , PO_4^{3-} , acetate, citrate, ethanol, glycerol, lactate, malate, oxalate

250 mg/L: NO_2^-

20 mg/L: SO_3^{2-} , ascorbate

1 mg/L: H_2O_2

0,1 mg/L: peracetic acid





Disposal:

Information regarding disposal can be found in the safety data sheet. You can download the SDS from www.mn-net.com/SDS.

Storage:

Avoid exposing the strips to sunlight and moisture. Keep container cool and dry (storage temperature 4–30 °C).

If correctly stored, the test strips may be used until the use-by-date printed on the packaging.



CTL SCIENTIFIC SUPPLY CORP. 1016-3 Grand Boulevard, Deer Park, NY 11729

Tel: 631-242-4249

Web: www.ctlscientific.com

Manufacturer: Macherey-Nagel GmbH & CO. KG

Rev: 2025-03