

REF 918128

Test 1-28

07.21

NANOCOLOR® iron LR**Method:**

Photometric determination of iron content by the triazine method.

Rectangular cuvette	50 mm	10 mm
Measuring range (mg/L Fe)	0.005–0.500	0.05–5.00
Measuring wavelength (HW = 5–12 nm)	563 nm / 540 nm	
Reaction time	3 min (180 s)	5 min (300 s)
Reaction temperature	20–25 °C	

Contents of reagent set

300 mL iron LR R1

20 g iron LR R2

1 measuring spoon 70 mm

Preliminary tests:

If there is uncertainty regarding the level of the concentration in the sample to be tested, a preliminary test with QUANTOFIX® Total Iron 100 (2–100 mg/L Fe, REF 91344) rapidly provides this information. This allows the dilution required for the determination to be calculated and prepared directly.

Interferences:

To test for the absence of interfering complexing agents we recommend a preliminary test with NANOCOLOR® organic complexing agent 10 (REF 985052).

The following will not interfere: ≤ 0.5 mg/L Co^{2+} ; ≤ 5 mg/L Al^{3+} , Mn^{2+} , MoO_4^{2-} ; ≤ 10 mg/L Cu^{2+} , Hg^{2+} , Ni^{2+} , Pb^{2+} , CN^- ; ≤ 50 mg/L $\text{Cr}_2\text{O}_7^{2-}$, NO_2^- ; ≤ 500 mg/L Ca^{2+} , Zn^{2+} ; ≤ 750 mg/L PO_4^{2-} ; ≤ 1000 mg/L Cd^{2+} , Mg^{2+} , NH_4^+ , SiO_3^{2-} ; < 1 % nonionic surfactants, cationic surfactants; ≤ 1 % anionic surfactants; ≤ 5 % sodium acetate; ≤ 20 % NaNO_3 .

* add 500 mg thiourea per 20 mL sample

The method is also suitable for the analysis of seawater. A prolonged reaction time of 5 min may be necessary.

Procedure:

Required accessories: Volumetric flask 25 mL (REF 91661), 10 mm and 50 mm rectangular cuvettes (REF 91933 and 91935), piston pipettes with tips

Sample	Blanc value ^[1]
In a 25-mL volumetric flask: Place 20 mL sample solution (the pH of the sample must be between pH = 1 and pH = 7), add 3 mL R1, mix, add 1 measuring spoon of R2	In a 25-mL volumetric flask, place: 20 mL distilled water, add 3 mL R1, mix, add 1 measuring spoon of R2
Fill up sample and blank value with distilled water to 25 mL and mix. After 3 min ^[2] pour into the cuvettes, clean the outside of the cuvettes and measure.	

^[1] In the iron determination in the 10-mm cuvette, the sample solution can also be used without addition of reagent as a blank value.

^[2] For the iron determination in the 10-mm cuvette a reaction time of 5 min is recommended.

Determination of iron(II) ions:

Procedure as described above, but without reagent ion R2.

Measurement:

For MACHEREY-NAGEL photometers see manual, test 1-28. Photometers of other manufacturers: Check the factor for each type of device by measuring standard solutions.

Analytical quality assurance:

NANOCONTROL Multistandard Metals 1 (REF 925015) or Multistandard Drinking Water (REF 925018)

Reduced analytical procedures:

To increase the number of determinations, 10 mL can be prepared in volumetric flasks: 8 mL sample solution + 1.2 mL R1 + 1 measuring spoon of R2. The evaluation is then performed in a semi-micro cuvette (REF 91950) .

Disposal:

Information regarding disposal can be found in the safety data sheet. You can download the SDS from www.mn-net.com/SDS.

CTL SCIENTIFIC SUPPLY CORP. 1016-3 Grand Boulevard, Deer Park, NY 11729

Tel: 631-242-4249

Web: www.ctlscientific.com

Manufacturer: Macherey-Nagel GmbH & Co. KG

REV: 2025-09