

pH-Fix 3.1–8.3

Test strips for measuring the pH values in the range of pH 3.1–8.3



Please refer to instructions for use

REF 92135

Pack content:

100 test strips

1 reference color scale pH 3.1–8.3

Application:

The pH-Fix test strips pH 3.1–8.3 are semi quantitative indicator test strips for many applications.

For quality control use commercially available pH buffer solution.

Instruction:

1. Dip the test strip with all its fields into the sample solution for 1 sec.
2. Shake off excess liquid.
3. Compare the colors with the color chart and read off the pH value which matches closest.



Remark:

In weakly buffered solutions, keep the test strip immersed for a longer period of time until color consistency of all test fields.

Storage:

Avoid exposing the strips to sunlight and/or moisture. Keep container cool and dry (store between 4–30 °C).

Disposal:

Used test strips may be disposed of as domestic waste.



Please read instructions for use!



Permitted storage temperature range



Use by



Batch identification



Item number



Manufacturer



Do not reuse