

REF 985007

Test 0-07 03.23

**NANOCOLOR® AOX 3**

Adsorbable organically bound halogens

**Method:**

AOX determination is carried out in 3 steps:

1. Solid phase extraction with **NANOSORB** for AOX
2. Decomposition of the concentrated adsorber medium
3. Determination as chloride with reagent set **NANOCOLOR® AOX 3**

Measuring range:	0.1 – 3.0 mg/L AOX	0.01 – 0.30 mg/L AOX
Method:	(0)071	(0)072
Wavelength (HW = 5 – 12 nm):	470 nm	
Reaction time:	3 min (180 s)	
Reaction temperature:	20 – 25 °C	

**Contents of the reagent set:**20 **NANOSORB** cartridges1 **preparation box** containing

- 2 × 100 mL rinsing solution concentrate for preparation of AOX 3 R1 (fill up each to 1 L with dist. water)
- 1 tube **NANOFIX AOX 3 R2**
- 1 × 105 mL AOX 3 R3
- 1 × 75 mL AOX 3 R4

20 **reaction tubes** 16 mm OD1 **detection box** containing

- 20 test tubes AOX 3
- 2 test tubes Chloride R2
- 1 test tube with blank value "NULL"

**Hazard warning:**Information regarding safety can be found on the box' label and in the safety data sheet. You can download the SDS from [www.mn-net.com/SDS](http://www.mn-net.com/SDS).**Interferences:**COD-contaminated test solution (> 50 mg/L COD) must be diluted beforehand. Alternatively, the **NANOCOLOR®** Supplement kit for AOX (up to 1000 mg/L COD, REF 918072) can be used.

When using 200 mL of rinsing solution, this method is also suitable for analyzing sea water.

**Procedure:**

Requisite accessories: Starter kit for AOX (REF 916111), piston pipette with tips; optional: supplement kit (REF 918072), pump set for AOX (REF 916115)

For exact measurements in the low range, it is recommendable using a real blank value:

**1a. Manual extraction**Connect a **NANOSORB** cartridge to the syringe 50 mL with the aid of an adaptor. Pour

**100 mL** test sample (*the pH value of the sample must be between pH 3 and 5*) into a glass beaker 150 mL, dip the **NANOSORB** cartridge into the test sample and lift the syringe plunger up and down 20 times to adsorb the organically bound halogens from the sample (accessories: stand with clamp and boss).

After extraction disconnect **NANOSORB** cartridge from the adaptor and syringe. Rinse the **NANOSORB** cartridge slowly in 4 – 5 portions with a total of

**100 mL R1** rinsing solution in order to remove inorganic chloride. Connect the syringe to the cartridge once more and blow out any excess of water from the **NANOSORB** adsorber with 2 strong draughts of air.


**1b. Extraction using the pump set**

Close valve of the flask. Pour

**100 mL** test sample (*the pH value of the sample must be between pH 3 and 5*) or **1000 mL** test sample (*the pH value of the sample must be between pH 3 and 5*) for the sensitive range into the flask and connect a **NANOSORB** cartridge to the flask using the adaptor. Open valve and start pumping for 20 min to adsorb the organically bound halogens from the sample.

After extraction disconnect the **NANOSORB** cartridge from the adaptor and flask. Rinse the **NANOSORB** cartridge in 4 – 5 portions with a total of

**100 mL R1** rinsing solution in order to remove inorganic chloride. Connect the syringe to the cartridge using the adaptor and blow out any excess of water from the **NANOSORB** adsorber with 2 strong draughts of air.

**2a. Decomposition if COD content is low, without supplement kit, using a heating block**

Add into a reaction tube 16 mm OD

1 **NANOFIX R2** and5 **mL R3** close and mix.

Open and insert the **NANOSORB** to this solution with help of a funnel, then press it down to the bottom of the tube with tweezers. Close the tube, place it into the heating block and heat at 120 °C for 30 min. Remove tube from heating block, shake gently and leave it to cool. Open tube, add

3.5 **mL R4**, close and mix.**2b. Decomposition if COD content is low, without supplement kit, using a microwave**

Add to the decomposition vessel

1 **NANOFIX R2** and5 **mL R3**, close and mix.

Open and add the **NANOSORB** to this solution using tweezers. Add a glass rod to the vessel to prevent the **NANOSORB** from swimming on the surface. Close the decomposition vessel. Place it on the outer edge of the microwave revolving plate and heat 23 s at 900 VA or 28 s at 750 VA (*always use the highest power rating of your microwave oven*).

Remove vessel from microwave and let cool for about 10 min. Turn the pressure vessel upside down once and open it with caution. Add

3.5 **mL R4**, close and mix.**2c. Decomposition if COD content is high, with supplement kit, using a heating block**

Add into a reaction tube 16 mm OD

1 **NANOFIX R2**,1 **black spoon R5** and5 **mL R3**, close and mix.

Open and insert the **NANOSORB** to this solution with help of a funnel, then press it down to the bottom of the tube using tweezers. Close the tube, place it into the heating block and heat at 120 °C for 30 min. Remove tube from heating block, shake gently and leave it to cool. Open tube, add

3.5 **mL R4** and

1 **orange spoon R6** (*the solution becomes turbid*), close and mix. Filter the solution with membrane or folded filters.

**2d. Decomposition if COD content is high, with supplement kit, using a microwave**

Add to the decomposition vessel

1 **NANOFIX R2**,1 **black spoon R5** and5 **mL R3**, close and mix.

Open and add the **NANOSORB** to this solution using tweezers. Add a glass rod to the vessel to prevent the **NANOSORB** from swimming on the surface. Close the decomposition vessel. Place it on the outer edge of the microwave revolving plate and heat 23 s at 900 VA or 28 s at 750 VA (*always use the highest power rating of your microwave oven*).

Remove the vessel from the microwave and let it cool for about 10 min. Turn the pressure vessel upside down once and open it with caution. Add

3.5 **mL R4** and

1 **orange spoon R6** (*the solution becomes turbid*), close and mix. Filter the solution with membrane or folded filters.

**3. Determination of AOX**

Open test tube AOX and add

4.0 **mL** decomposition solution (*let particles of adsorbent deposit or use membrane filters*). Add1.0 **mL Chloride R2**, close and mix.

Clean outside of test tube and measure after 3 min.

Adjust photometer to zero by using blank value "NULL".

**Measurement:**

For using MACHEREY-NAGEL photometers see manual, test 0-07.

**Photometers of other manufacturers:**

For other photometers check whether measurement of round glass tubes is possible. Verify calibration curve for each type of instrument by measuring standard solutions.

**Analytical quality control:****NANOCONTROL AOX 3** (REF 92507)

The measurement of a blank value and a standard is recommended before every measuring series as quality control measure.

**Reference:**

German Standard Methods for the Examination of water, waste water and sludge (DIN EN 9562)

CTL SCIENTIFIC SUPPLY CORP. 1016-3 Grand Boulevard, Deer Park, NY 11729

Tel: 631-242-4249

Web: [www.ctlsupply.com](http://www.ctlsupply.com)

Manufacturer: Macherey-Nagel GmbH &amp; CO. KG

Rev: 2025-03