

# Nitratesmo

## Test Paper for the rapid determination of Nitrate and Nitrite

### Colour reaction:

The white test paper turns red-orange in the presence of nitrate; lemon-yellow in the presence of nitrite.

### Method of Application:

#### a) Detection of Nitrate

Dip test paper briefly into test solution and apply excess liquid to a strip of filter paper. Dip the moistened test paper very briefly into conc. sulphuric acid and spread on a glass plate placed on white paper. In the presence of nitrate only, the paper turns red. If, in addition to nitrate, nitrite is also present (evidenced by yellowish-red colour), the interfering nitrite can be destroyed by using amido-sulphuric acid.

**Limit of sensitivity:** 10 mg/l  $\text{NO}_3^-$ .

#### b) Detection of Nitrite

Dip test paper briefly into test solution and apply excess liquid to a strip of filter paper. Dip the moistened test paper very briefly into 5N hydrochloric acid (17%) and spread on a glass plate placed on white paper. The presence of nitrite is indicated by a yellow colour, which, in the case of minute quantities, appears only after 3-4 minutes.

**Limit of sensitivity:** 5 mg/l  $\text{NO}_2^-$ .

### Note:

For the determination of nitrate or nitrite in concentrations approaching the limit of sensitivity, a control test is recommended. Do not dip test paper strip too frequently into the same sulphuric or hydrochloric acid.

### Interferences:

Chlorate, bromate, iodate, hypochlorite, vanadate and iodide interfere with the Nitratesmo reaction.

#### Product data and ordering information

<b>Type</b>	<b>Qualitative test paper</b>
<b>Presentation</b>	<b>Reel of 5 m length and 10 mm width</b>
<b>Virage</b>	<b>Nitrate: white to red</b> <b>(dip test paper into sample and afterwards into sulfuric acid 96%)</b> <b>Nitrite: white to yellow</b> <b>(dip test paper into sample and afterwards into Hydrochloric acid 5 mol/l)</b>
<b>Limit of sensitivity</b>	<b>10 mg/l (ppm) Nitrate (<math>\text{NO}_3^-</math>),</b> <b>5 mg/l (ppm) Nitrite (<math>\text{NO}_2^-</math>)</b>
<b>REF</b>	<b>90611</b>

CTL SCIENTIFIC SUPPLY CORP. 1016-3 Grand Boulevard Deer Park, NY 11729

Tel: 631-242-4249

web: [www.ctlscientific.com](http://www.ctlscientific.com)

Manufacturer: Macherey-Nagel GmbH & CO. KG