

# DIPYRIDYL Paper

for the rapid determination of ferrous ions

## Colour reaction:

In the presence of Fe(II)-ions a red spot appears on white background.

## Method of application:

Apply a drop of the weak mineral acid test solution (> pH 2) to the test paper. In the presence of Fe(II)-ions a red spot appears, or in the presence of small quantities of Fe(II)-ions a red ring.

**Limit of sensitivity:** 2 mg/l Fe(II)

## Interferences:

The paper is specific for the determination of Fe(II)-ions. A few other metals result in pale discolourations of the test paper, but these are covered by the colour reaction of the ferrous ions. Fe(III)-ions result in a pale yellow reaction. In the presence of substantial quantities of Fe(III)-ions and small quantities of Fe(II)-ions, the Fe(III)-ion can be masked by the addition of small quantities of NaF.

### Product data and ordering information

Type	Qualitative test paper
Presentation	Pack of 200 strips 20x70 mm
Color reaction	white to red
Limit of sensitivity	2 mg/l (ppm) Fe <sup>2+</sup>
REF	90725

CTL SCIENTIFIC SUPPLY CORP. 1016-3 Grand Boulevard Deer Park, NY 11729

Tel: 631-242-4249

web: [www.ctlscientific.com](http://www.ctlscientific.com)

Manufacturer: Macherey-Nagel GmbH & CO. KG