

Chlorine Test

Chlorine Test is a test paper for the semi-quantitative determination of total chlorine and free chlorine in the concentration range of 10–200 mg/L.

Color reaction:

The white test paper changes to blue-violet if chlorine is present.

Instructions for use:

Briefly immerse the test strip in the test solution and shake off excess liquid. **Immediately** compare reaction color with the color scale.

Note:

If the reaction color disappears after immersion, chlorine concentrations higher than 200 mg/L are present. In this case, the result should be discarded.

Storage conditions:

Protect test paper from sunlight and moisture. Store package in a cool and dry location (storage temperature 4-30 °C).

CTL SCIENTIFIC SUPPLY CORP. 1016-3 Grand Boulevard, Deer Park, NY 11729

Tel: **631-242-4249**

Web: www.ctlscientific.com

Manufacturer: **Macherey-Nagel GmbH & CO. KG**

Rev: 2025-03